























QUICK REFERENCE GUIDE

1>> Open Futudent Recording Software to activate live streaming	
	The active camera will livestream without recording . Position the camera over the treatment area.
	Place the wireless foot pedal in position for the person who will record the images.
2>> Record video and set bookmarks (2 options)	
	On-Screen - click the recording icon to start or stop recording. To Add Bookmark  during recording press the bookmark icon briefly.
	Foot pedal - press the pedal for 3 seconds to start or stop recording. To Add Bookmark  during recording press the pedal briefly.
3>> Take still images (2 options)	
	On-screen - click the camera icon
	Foot pedal – briefly press the foot pedal
4>> Open Gallery	
	Click the arrows on the left side of the screen to open and close the Gallery .
 Videos  Photos	Select Videos or Photos in the Gallery to see case documentation.
 Gallery  Camera	Go to the Gallery or Camera livestream by pressing the buttons at the bottom-left of the screen.
5>> Add annotations and more bookmarks	
	Open a video by pressing the play arrow on the video thumbnail in the gallery.
 	Move the video time elapsed bar to the desired location and press the Add Bookmark  icon. The bookmark appears as a line.
	Click the arrows on the right side of the video to open and close the Bookmarks .
	Double click a bookmark thumbnail to open the annotation screen to add notes.
6>> Save the recording or images on the clinic server or computer	
	Still images can be copied and pasted in to other imaging software or patient records with import functionality.
7>> Share with colleagues or patients using the Cloud or MyDentalBook	
	Send to patient. Attach treatment plans and home-care instructions.
	Send to colleague. Your colleague will receive a link to follow.
	Upload to cloud. Save your cases for future consultation.