

AXIS™ LED LIGHT



Axis™ LED Light

Activ™ Neck Battery

Precision Dentistry

A light attached to your loupe leverages coaxial alignment to remove shadows that overhead light causes

Visibility Without Patient Discomfort

Loupe lights cast a spot wide enough to fill your loupe image without shining an uncomfortable light in your patients eyes

Increase Resolution

Illumination adds yet another level of resolution, meaning you are able to make the best clinical decisions given patient need

Advanced Comfort and Visualization

Enova developed the **Axis™ LED Clip-On Light** to be a lightweight, ergonomic illumination solution compatible with all Enova loupes. An accompanying neckband-style battery powers 30,000 lux illumination, variable light intensity, maximum magnification field of view, and 5+ hours of working time. When used in conjunction with lightweight Vizix® loupes, clinicians benefit from superior illumination, visualization, and comfort, allowing them to perform with ease and confidence.

ERGONOMIC NECK BATTERY

Cut out the hassle of struggling with traditional tethered battery packs! The neck-worn, capacitive touch Activ battery can be easily activated and adjusted through disposable gowns and without removing gloves.

COLOR PERSONALIZATION

Choose from five different distinct color caps to personalize your light. Available in silver, black, blue, green and purple. Activ™ Battery caps also come in black, green, blue and purple to complete the look.

OPTIONAL HD GLASS COMPOSITE FILTER

Removes the light wavelength that tends to speed up the curing of composite materials specific to the dental field. When not needed, simply flip the filter up for pure LED light.

TECHNICAL DATA

Light Source: High Performance LED
 Lamp Illumination: 30,000 lux
 Color Temperature: 6,100k
 Illumination Spot Size: 4 inches, fixed
 Optional: HD glass composite filter
 Construction: Aircraft Aluminum
 Weight: 7 grams
 LED Warranty: 3 years
 Battery: Rechargeable Lithium-Ion (charger included)
 Battery Warranty: 1 year

| | | | | |
|--------------------------------|----------------------------|---|--------------|---|
| Battery Pack Options | Activ Neck Battery* |  | PP1.5 |  |
| Battery Run-Time | 4-8 hours | | 10-20 hours | |
| Battery Pack Dimensions | 21" long, 1"x1" | | 4"x2.5"x1" | |
| Battery Pack Weight | 5oz (149g) | | 7oz (200g) | |

*Recommended by Enova

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of LED clip-on lights and surgical headlamps so you can select the specifications to fit your needs. Whether you need bright, penetrating light for deep surgical cavities or cordless freedom for your active work, Enova has you covered. Visit www.goenova.com to learn more.



Qubit Cool (XLT-105F)
 105,000 lux
 Tethered wearable battery
 Run-time: 8+ hours



Axis Light™
 30,000 lux
 Wearable neck battery/
 traditional battery pack
 10+ hour run-time



Atom™ Light
 23,000 lux
 Wireless Battery
 2+ hour run-time

To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
 6305 Sandburg Road STE 100
 Minneapolis, MN 55427

