

proCAM XS BY FUTUDENT



Versatile mount is compatible with most loupes, headlights, and overhead lights

Optional customized light adapter

Adjust focus to your exact working distance



Lightweight: At just 14 grams, you'll barely notice a difference when mounted to your loupe or headlight

Durable Aluminum Housing

LIGHTS, CAMERA, DOCUMENTATION

Modern healthcare presents an increasing need for safe, non-invasive and efficient video and photo capture for clinical documentation, lecturing, and education. Weighing just 14 grams, the **Futudent proCam XS** is a miniature camera with giant 4K features and easy mounting to loupes, headlights and overhead lights.

FAST, SIMPLE EDITING

Proprietary Futudent software captures and instantly makes your video and photos available so you can focus on your patient, not your camera.

HANDS-FREE OPERATION

Start and stop video recording and capture photos without lifting a finger. The proCam XS is powered by a footpedal so you don't have to interrupt treatment to document it.

TREATMENT ACCEPTANCE

Increase your ROI and treatment acceptance rates by being able to instantly help patients visualize their case and comprehend your treatment plan.

TECHNICAL DATA

Resolution: 4K30/1080p30/720p30/8Mpix stills

Weight: 14 grams (0,49 ounces)

Dimensions: 1.18" x 0.79" x 0.79" (30 × 20 × 20 mm)

Optics: Exchangeable UHD lenses with manual focus (8/16/25/35 mm)

Construction: Full Aluminum Body

Connectivity: USB 2 with detachable cable (USB micro connector)

Video Editing Software: Futudent software; 5 licenses included/unit


Camera Warranty: 2 years

Compliance: MDR 2017/745, Classification 1

Mounting Options	
Loupe Mount	Universal and Specific Loupe Adapters
LED Headlight	Headband and Headlight Adapters
Free Standing	Overhead Light, Gooseneck/Flexible Arm Adapters

PAIRS WITH YOUR FAVORITE ENOVA PRODUCTS

The latest cameras from Futudent are so lightweight, you'll hardly notice the difference when you mount the camera on your loupe or headlight.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road STE 100
Minneapolis, MN 55427

F20474, Rev B
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

