

SMARTCAM XS BY FUTUDENT



LIGHTS, CAMERA, DOCUMENTATION

On the leading edge of healthcare documentation technology, the **Futudent smartCam** is the **world's first auto-focusing medical camera**. This 4k miniature camera is a great solution for high level documentation, communication, and live streaming for communication and educational purposes.

FAST, SIMPLE EDITING

Proprietary Futudent software captures and instantly makes your video and photos available so you can focus on your patient, not your camera.

HANDS-FREE OPERATION

Start and stop video recording and capture photos without lifting a finger. The smartCam is powered by a footpedal so you don't have to interrupt treatment to document it.

TREATMENT ACCEPTANCE

Increase your ROI and treatment acceptance rates by being able to instantly help patients visualize their case and comprehend your treatment plan.

TECHNICAL DATA

Resolution: 8.3 MP (4k: 3840×2160)

Weight: 50 grams (1.76 ounces)

Dimensions: 1.57" x 0.87" x 1.06" (40 × 22 × 27 mm)

Optics: 25mm Novocam UHD Lens

Construction: Full Aluminum Body

Connectivity: Micro USB female on camera

Video Editing Software: Futudent software; 5 licenses included/unit

Camera Warranty: 2 years

Compliance: MDR 2017/745, Classification 1

Mounting Options	
Loupe Mount	Universal and Specific Loupe Adapters
LED Headlight	Headband and Headlight Adapters
Free Standing	Overhead Light, Gooseneck/Flexible Arm Adapters

PAIRS WITH YOUR FAVORITE ENOVA PRODUCTS

The latest cameras from Futudent are so lightweight, you'll hardly notice the difference when you mount the camera on your loupe or headlight.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road STE 100
Minneapolis, MN 55427

F20504, Rev A
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

