

QUASAR™ COOL LED HEADLIGHT



Adjustable light angle
for ergonomic wear

Durable Aluminum Housing

Adjustable spot

Adjust fit from over
and behind the head
for comfort and no-
slip wear

All headband surfaces
treated with Aegis®

THE MOST VERSATILE LED HEADLIGHT

The illumination spot of the **Quasar Cool LED Headlight** features clear, cool white light that is uniform and consistent. This headlight offers the adjustability that makes it our most versatile and popular surgical headlight.

150,000 - 300,000 LUX

Lux* ranging from 150,000 - 300,000 — depending on desired working distance — of pure white homogenous light for extreme surgical field illumination.

COMFORTABLE AND ERGONOMIC

Lightweight construction, ergonomic and tangle-free battery connectivity, and fit-adjust features make this headlight comfortable for long-term wear.

INFECTION PREVENTION

Enova headlights do not require fans, all headbands are treated with antimicrobial coating, and easy-to-clean light housings simplify infection prevention and control.

*Lux is variable based on preferred working distance.

TECHNICAL DATA

Light Source: LED refractive technology

Lamp Illumination: 150,000 (at 14") - 300,000 (at 10") lux

Color Temperature: 6,100K

Illumination Spot Size: Adjustable from 2" to 5" (5.08cm to 12.7cm)

Headlight Construction: All Aluminum

Headlight Weight: 6oz (170g)

Headband Construction: Nylon Copolymer, Replaceable Comfort Pads

LED Warranty: 5 years

Battery: Rechargeable Lithium Ion (base station charger included)

Battery Warranty: 1 year

Battery Pack Options	PP4	PP2	FREESTYLE
Battery Run-Time	10-24 hours	5-10 hours	4-8 hours
Battery Pack Dimensions	3"x4"x1"	4"x2.25"x1"	5.5"x2.375"x.0875"
Battery Pack Weight	9 oz (255g)	7oz (200g)	6oz (170g)

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of the brightest, lightest surgical headlights so you can select the exact specifications to fit your needs. Try an Enova headlight with our free evaluation program! Visit www.goenova.com to learn more.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road STE 100
Minneapolis, MN 55427

F20430, Rev B
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Microban® is a registered trademark of Aegis, Inc.



QUASAR™ WARM LED HEADLIGHT



Adjustable light angle
for ergonomic wear

Durable Aluminum Housing

Adjustable Spot

Adjust fit from over
and behind the
head for comfort
and no-slip wear

All headband surfaces
treated with Aegis®

SEE BETTER. PERFORM BETTER.

With color-rendering index (CRI) of 95, the **Quasar Warm LED Headlight** offers the truest color rendition on the market. While providing slightly less brightness than the Enova Q-Cool series, these headlights illuminate in near true-color for the greatest tissue differentiation.

PURE LIGHT TECHNOLOGY

Produces best-in-class color rendition for greater tissue distinction while reducing eye strain with lux* of 115,000 - 230,000.

COMFORTABLE AND ERGONOMIC

Designed for comfortable, long-term wear, the Quasar Warm headlight has several adjustable features to help you shape the fit and angle to your exact needs.

INFECTION PREVENTION

Enova headlights do not require fans, all headbands are treated with antimicrobial coating, and easy-to-clean light housings simplify infection prevention and control.

*Lux is variable based on preferred working distance.

TECHNICAL DATA

Light Source: LED refractive technology

Lamp Illumination: 115,000 (at 14") - 230,000 (at 10") lux

Color Temperature: 4,500K

Illumination Spot Size: Adjustable from 2" to 5" (5.08cm to 12.7cm)

Headlight Construction: All Aluminum

Headlight Weight: 6oz (170g)

Headband Construction: Nylon Copolymer, Replaceable Comfort Pads

LED Warranty: 5 years


Battery: Rechargeable Lithium Ion (base station charger included)

Battery Warranty: 1 year

Battery Pack Options	PP4	PP2	FREESTYLE
Battery Run-Time	10-24 hours	5-10 hours	4-8 hours
Battery Pack Dimensions	3"x4"x1"	4"x2.25"x1"	5.5"x2.375"x.0875"
Battery Pack Weight	9 oz (255g)	7oz (200g)	6oz (170g)

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of the brightest, lightest surgical headlamps so you can select the exact specifications to fit your needs. Try an Enova headlight with our free evaluation program! Visit www.goenova.com to learn more.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandbrug Road STE 100
Minneapolis, MN 55427

F20429, Rev B
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Microban® is a registered trademark of Aegis, Inc.

