

QUBIT™ COOL LED HEADLIGHT



Adjustable light angle
for ergonomic wear

Adjust fit from over
and behind the
head for comfort
and no-slip wear

Durable Aluminum Housing

Fixed Spot

**NEW! Optional attachable Beam Adjustment Filter
doubles the otherwise fixed spotlight size**

All Enova LED headlights are
manufactured in the USA

THE MOST ECONOMICAL LED HEADLIGHT

The illumination spot of the **Qubit Cool LED Headlight** features clear, cool, white light that is uniform and consistent. This headlight offers a lightweight and economical way to bring LED light into any practice. Attach the the Beam Adjustment Filter to make the headlight semi-adjustable.

105,000 - 210,000 LUX

Lux ranging from 105,000 - 210,000 — depending on desired working distance — of pure white homogenous light for extreme surgical field illumination while feather light.

COMFORTABLE AND ERGONOMIC

Lightweight construction, ergonomic and tangle-free battery connectivity, and fit-adjust features make this headlight comfortable for long-term wear.

BEAM ADJUSTMENT

The simple flip filter allows the Qubit to have a semi-adjustable light that doubles the spot size. For smaller, more focused light, simply adjust the filter off the lens.

TECHNICAL DATA

Light Source: LED refractive technology

Lamp Illumination: 105,000 (at 14") - 210,000 (at 10") lux

Color Temperature: 6,100K

Illumination Spot Size: Fixed spotlight at 2" (5.08cm), Optional Beam Adjustment Filter

Headlight Construction: All Aluminum

Headlight Weight: 7 oz/198g

Headband Construction: Nylon Copolymer, Replaceable Comfort Pads

LED Warranty: 5 years

Battery: Rechargeable Lithium Ion (base station charger included)

Battery Warranty: 1 year

Battery Pack Options	PP4	PP2	FREESTYLE
Battery Run-Time	10-24 hours	5-10 hours	4-8 hours
Battery Pack Dimensions	3"x4"x1"	4"x2.25"x1"	5.5"x2.375"x.0875"
Battery Pack Weight	9 oz (255g)	7oz (200g)	6oz (170g)

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of the brightest, lightest surgical headlamps so you can select the exact specifications to fit your needs. Try an Enova headlight with our free evaluation program! Visit www.goenova.com to learn more.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road STE 100
Minneapolis, MN 55427

F20437, Rev B
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Microban® is a registered trademark of Aegis, Inc.



QUBIT™ WARM LED HEADLIGHT



Adjustable light angle for ergonomic wear

Durable Aluminum Housing

Fixed Spot

Adjust fit from over and behind the head for comfort and no-slip wear

NEW! Optional attachable Beam Adjustment Filter doubles the otherwise fixed spotlight size

All Enova LED headlamps are manufactured in the USA

THE MOST ECONOMICAL LED HEADLIGHT

With color-rendering index (CRI) of 95, the **Qubit Warm LED Headlight** features a warm light that offers the truest color rendition on the market. This headlight offers a lightweight and economical way to bring LED light into any practice. Attach the Beam Adjustment Filter to make the headlight semi-adjustable.

80,000 - 160,000 LUX

Lux ranging from 80,000 - 160,000 — depending on desired working distance — of pure homogenous light for true tissue color rendition.

COMFORTABLE AND ERGONOMIC

Lightweight construction, ergonomic and tangle-free battery connectivity, and fit-adjust features make this headlight comfortable for long-term wear.

BEAM ADJUSTMENT

The simple flip filter allows the Qubit to have a semi-adjustable light that doubles the spot size. For smaller, more focused light, simply adjust the filter off the lens.

TECHNICAL DATA

Light Source: LED refractive technology

Lamp Illumination: 80,000 - 160,000+ lux

Color Temperature: 4,500K

Illumination Spot Size: Fixed spotlight at 2" (5.08cm), Optional Beam Adjustment Filter

Headlight Construction: All Aluminum

Headlight Weight: 7 oz/198g

Headband Construction: Nylon Copolymer, Replaceable Comfort Pads

LED Warranty: 5 years

Battery: Rechargeable Lithium Ion (charger included)

Battery Warranty: 1 year

Battery Pack Options	PP4	PP2	FREESTYLE
Battery Run-Time	10-24 hours	5-10 hours	4-8 hours
Battery Pack Dimensions	3"x4"x1"	4"x2.25"x1"	5.5"x2.375"x.0875"
Battery Pack Weight	9 oz (255g)	7oz (200g)	6oz (170g)

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of the brightest, lightest surgical headlamps so you can select the exact specifications to fit your needs. Try an Enova headlight with our free evaluation program! Visit www.goenova.com to learn more.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandberg Road STE 100
Minneapolis, MN 55427

F20438, Rev B
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Microban® is a registered trademark of Aegis, Inc.

